



Enhance Feeding Pump

General Features:

- ▶ Pump supports all kind of nutrition and liquid foods.
- ▶ 25 ml/hr to 600 ml/hr flow rate.
- ▶ Target volume selection available from 1 ml to 1500 ml.
- ▶ Automatically stops the pump reaching the target volume.
- ▶ Accuracy: +/- 10%.
- ▶ Last position memory retention available.
- ▶ Flashing display with alarm indication.
- ▶ Display of calorie value.

Technical specification:

- ▶ Power supply DC 19V, 3-5A.
- ▶ Battery 14.8v rechargeable sealed Nickel Metal Hydride.
- ▶ Display LCD and seven segment display
- ▶ Flow rate range: 25 ml/hr to 600 ml/hr
- ▶ Target volume range: 1 ml to 1500 ml.
- ▶ Elapsed time: Display of Elapsed time from 00.00 to 24.00 hrs.
- ▶ Dimension: 140(L) × 145(B) × 205(H)mm.
- ▶ Weight: 2.3 Kg with Pole Clamp.

Modes:

- ▶ Prime mode.
- ▶ Normal mode.
- ▶ Flush mode: 1000 ml/hr

Alarm indication:

- ▶ Bag empty.
- ▶ Target volume reached.
- ▶ Prime mode completed.
- ▶ Flush mode completed.
- ▶ Key disconnected.
- ▶ Validity of key expired.
- ▶ Save parameters.



Note: Enhance Feeding Bag can only be used.

AKAS Medical

240/1, Periya Colony, Athipet, Ambattur, Chennai - 600 058
Cell: +91-98403 79116, Tel: +91-44-3253 3333, Fax: +91-44-2635 0030
enquiry@akasmedical.com www.akasmedical.com

*Specifications Subject to change without notice.

Your nearest dealer: